

Evaluation of Drug Mediated Changes in Action Potentials & Calcium Transients Recorded from Adult Human Stem Cell-derived Cor.4U[®] Cardiomyocytes

Ionic Transport Assays

Ken Gibson, Ph.D., Yimei Yue, Ping Li, Ph.D.
Ionic Transport Assays, Inc., Saint Louis, Missouri, USA

Ionic Transport Assays


Abstract

- The FDA CiPA initiative proposes the use of adult human induced pluripotent stem cell-derived cardiomyocytes for testing proarrhythmia. This work contrasts the actions of diverse compounds on action potentials and internal calcium transients recorded from Cor.4U[®] Cardiomyocytes.
- Prolongation of action potential duration (APD₉₀), early after depolarizations (EADs) and decreased spontaneous beating rate of calcium transients were observed after exposure to the potassium channel blockers, E-4031, dofetilide and sotalol.
- Reduction of APD₉₀, shortened duration and increased spontaneous beating rate of internal calcium transients were observed after exposure to the calcium channel blockers, nifedipine and verapamil.
- Prolongation of action potential duration (APD₉₀) was observed after exposure to the sodium channel blocker, flecainide.
- Cor.4U[®] Cardiomyocytes provide a robust model for cardiac safety pharmacology.


Methods

- Adult human heart cell generation and cell culture**
For this study, single vials containing ~4.5 x 10⁶ Cor.4U[®] cardiomyocytes (Axiogenesis AG) were thawed and plated into 6 well plates treated with 0.1% gelatin. This was defined as culture day 1 for the purpose of this study. Experiments were performed on Cor.4U[®] cardiomyocytes after 10 - 20 days of growth in the incubator (7% CO₂ and 37°C).
- Manual Current Clamp Action Potential Assay**
Perforated patch clamp recordings were performed with the external perfusate solution composed of: 150 mM NaCl, 5.4 mM KCl, 1.8 mM CaCl₂, 1.0 mM MgCl₂, 15 mM glucose, 15 mM HEPES, 1 mM Na-pyruvate (pH 7.4) and the internal, or pipette solution contained: 150 mM KCl, 5 mM NaCl, 2 mM CaCl₂, 5 mM MgCl₂, 5 mM EGTA, 10 mM HEPES (pH 7.2). Cor.4U[®] cardiomyocytes were paced at 1 Hz and the increasing concentrations of test compounds were administered at four minute intervals. If significant increases in action potential duration were observed, the external pacing was stopped and spontaneous action potential was observed for EADs.
- HTS Internal Calcium Transient Assay**
On culture day 5-7, Cor.4U[®] cardiomyocytes were re-suspended with trypsin and re-plated at a desired density (>10,000) in a black clear-bottom 96 well plate which were pre-coated with 0.1% gelatin. After another 10 - 15 days of growth, EarlyTox[™] Cardiotoxicity Kit (Molecular Devices) was used to perform calcium transient assay. A baseline measurement was recorded and reference test compounds were administered to each plate. The plates were maintained in the incubator for 30 min and then transferred to the SpectraMax i3 for fluorescence reading at 37°C.

Cor.4U[®] Cardiomyocyte Interaction of the Action Potential and Cardiac Contraction



Calcium Fluorescence* Indicator




- Peak Count & Frequency
- Peak Position (time) and Amplitude
- Peak Width (FWHM)
- Rise Time (10% to 90%)
- Decay Time (90% to 10%)


*SpectraMax i3, SoftMax Pro 6.4, Molecular Devices, Sunnyvale, CA

- An action potential is the interactive summary of transmembrane ion currents to produce electrical activation of the myocyte leading to an increase in cytoplasmic free Ca²⁺, the internal biochemical stimulation of myocyte contraction.


Potassium Channel Blockers Increase APD₉₀ in Ventricular-like Cor.4U[®] Cells




Nifedipine Blocks I_{Ca++}, Shortens APD & Increases Spontaneous Beating Rate



Verapamil Shortens, then Increases APD₉₀



Flecainide Blocks I_{Na+} & Prolongs APD₉₀



Summary of Observed Effects

Compound	Manual Electrophysiology: 1 Hz				Spontaneous Beating Frequency (BPM)
	APD ₆₀ (ms)	APD ₉₀ (ms)	Rise-rate (v/s)	RMP (mv)	
E-4031	Red	Red	Red	Yellow	Red
Dofetilide	Red	Red	Red	Yellow	Red
Terfenadine	Green	Green	Green	Yellow	Red
Ibutilide	Green	Green	Green	Yellow	Red
Moxifloxacin	Green	Green	Green	Yellow	Red
Cisapride	Green	Green	Green	Yellow	Red
Astemizole	Green	Green	Green	Yellow	Red
Sotalol	Green	Green	Green	Yellow	Red
Nifedipine	Green	Green	Green	Yellow	Green
Nitrendipine	Green	Green	Green	Yellow	Green
Verapamil	Green	Green	Green	Yellow	Green
Flecainide	Green	Green	Green	Yellow	Red
Mexiletine	Green	Green	Green	Yellow	Red
Quinidine	Green	Green	Green	Yellow	Green
Risperidone	Green	Green	Green	Yellow	Not Tested
Stausrosprine	Green	Green	Green	Yellow	Not Tested
Chlorpromazine	Green	Green	Green	Yellow	Not Tested
ATX-II	Green	Green	Green	Yellow	Not Tested

Decrease ≥ 50% Decrease > 20-49% No change Increase > 20-49% Increase ≥ 50%



Ionic Transport Assays, St. Louis, MO
www.ionictransport.com | 314-743-7922